



# FloorBed<sup>®</sup> 2

By Accora

## Specification Sheet



**Reduce the risk of injury in the event of a patient fall**

**Reduce the risk of caregiver injury**

**Enable independence and mobility**

# OVERVIEW

**Order code:**

**NSB-0-FL2-200**

# SPECIFICATION

Overall dimensions:	925 mm x 2255 mm (36.4 x 88.8 in)
Mattress size (standard configuration):	900 mm x 2000 mm (35.4 x 78.7 in)
Mattress size (extended configuration):	900 mm x 2200 mm (35.4 x 86.6 in)
Safe working load:	185 kg (29 st)
Patient weight:	150 kg (24 st)
Lift range (excl. mattress):	71 mm - 800 mm (2.8 in – 31.5 in)
Profiling functions:	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg

# TECHNICAL INFORMATION

<b>Head and foot bumpers</b> BMHNE-0-FL2-100	<ul style="list-style-type: none"> <li>Consists of both foot and head bumper.</li> <li>Attachable to the frame by means of nylon straps and plastic buckle clips.</li> </ul>
<b>Standard bed lever</b> STLEV-0-FL1-100	<ul style="list-style-type: none"> <li>Can be placed in several different positions on the bed frame.</li> <li>Provides support for bed mobility and transfers.</li> </ul>
<b>Rotating bed lever</b> RTLEV-0-FL1-000	<ul style="list-style-type: none"> <li>The additional length provides extra support with lying to sitting movement.</li> <li>When rotated, this provides additional support from sit to stand, and may enable to "pull" themselves up rather than push.</li> </ul>
<b>Integrated padded side rails – standard length</b> SDR-0-FL2-000	<ul style="list-style-type: none"> <li>Full length support.</li> <li>Easy lock and lower mechanism that locks securely in position.</li> <li>Side rails do not impede on transfers as they are below the mattress when lowered.</li> <li>2 in 1 bumper and side rail.</li> <li>Compatible with headboard/ footboard bumpers.</li> </ul>
<b>Integrated padded side rails – for use when FloorBed® is extended by 20 cm</b> SDREX-0-FL1-000	<ul style="list-style-type: none"> <li>Total length 220 cm.</li> </ul>
<b>FloorBed® safety sleeve</b> SLE-0-FL1-000	<ul style="list-style-type: none"> <li>Safety sleeve may reduce risk of entrapment by reducing gaps when raising the knee break and back rest functions (protecting the platform mechanism).</li> <li>Risk of entrapment is further reduced when connected to the footboard and headboard bumpers.</li> </ul>
<b>Support bar / Self assist pole</b> LIFOL-0-FL1-000	<ul style="list-style-type: none"> <li>Support for sitting up or changing position in bed.</li> </ul>
<b>High safety mat</b> SAFTH-0-FL1-000	<ul style="list-style-type: none"> <li>22 cm high safety mat - top of safety mat is same height as the top of the mattress when a 15 cm mattress is used, and FloorBed® is at its lowest position.</li> <li>Folds for storage.</li> </ul>
<b>High safety mat with integrating sliding sheets</b> SAFTHF-0-FL1-000	<ul style="list-style-type: none"> <li>22 cm high safety mat - top of safety mat is same height as the top of the mattress when a 15 cm mattress is used, and FloorBed® is at its lowest position.</li> <li>Integrated slide sheets for ease of transfer.</li> <li>Folds for storage.</li> </ul>
<b>Multicare 3000 Visco Mattress</b> VISM-0-FM1-000	<ul style="list-style-type: none"> <li>Layered soft memory foam technology for improved pressure reduction.</li> <li>Base layer of high density foam for support.</li> <li>Robust foam structure for longevity.</li> </ul>
<b>Mattress infill</b> MTINF-0-FL1-000	<ul style="list-style-type: none"> <li>Infills to fit between mattress and side rails.</li> </ul>